Technical specifications

Height	68 / 89 cm
Width	76 cm
Total length	185 cm
Base dimensions	121x54 cm
Cushion thickness	13 cm
Adjustment of the backrest angle	0°~+70°
Adjustment of the legrest angle	0°~-20°
Adjustable armrests	Yes, in height (±10 cm)
Facial hole	Yes
Maximum loading weight	200 kg
Upholstery	PU

Number of motors Hand controller Wheels Reference Voltage Frequency Consumption Net weight Gross weight Packaging size Assembly

4 Yes (standard) Optional 2249 / 2249 brown base 220 V /240 V 50/60 Hz 0,12 kWh 104 kg 124 kg 204x81x72 cm It is delivered assembled

Performances

By means of the standard hand controller or the optional foot controller, almost all the table movements can be adjusted.



The large height-stroke (68 cm to 89 cm) eases the accessibility of the patients to the table.



The legrest and backrest inclination allows a wide range of treatments with the patient face up.



The anti-Trendelenburg position helps the exit from the table to the patients with reduced mobility or elderly people.



The Trendelenburg position is commonly used for lymph draina-

ge to improve blood flow and fluid

retention on the legs.



The headrest is equipped with a facial hole to provide comfort to the patients when they are face down.



The armrests can be adjusted by a gaslift and they also help the professional during facial massages.



Position adjustments



If you press the first button, the backrest inclination can be adjusted. If you press the second button, the legrest inclination can be adjusted. If you press the third button, the inclination of the whole table can be adjusted, in the Trendelenburg position. If you press the fourth button, the inclination of the whole table can be adjusted, in the

Antitrendelenburg position. If you press the fifth button, the height of the table can be adjusted.